



NEO Apex Locator

3 Frequencies
World's Smallest Apex Locator

Accurate Technology for Endodontic Treatment



Our Polaroid Apex Locator Include Six Various Accessories.

Polaroid

Website: polaroidhealth.com

E-Mail: info@polaroidhealth.com

Phone: 1-844-789-5050



Small and Compact

Slick touch screen design, with a compact size allowing you to be mobile and ready when needed. No need for bulky hard to use machines.



Various Attachments

Attach the device on a stand, your apron, or a dental ruler that is included with your purchase. (Now comes with a finger clip to avoid having to move your head back and forth during treatment.)



Touch OLED Screen

Not just a display, but the focal point of the whole experience, improve treatment productivity with our bright and clear OLED screen display. (3 display programs available – For your conveniences.)



Calibration Made Simple

Calibrate the device by yourself with our calibration jig. Calibration ensures accuracy and eliminates the effect of changes in temperature, moisture etc. inside the canal even during treatment.



W.L.S Mode

(Working Length Select) Adjustable working length for ease of use. Select your convenience with working length setting, determine the root canal working length during root canal treatment.



Included with Unit

One Probe Cord, One Extension Cord, two File Hooks and Three Lip Hooks Included with every unit.



Auto On/Off

Device automatically turns off after 3 minutes if not in use. Conserve's power in the long run.



Accurate Measurement

Using 3 frequencies for better results and the latest signal IC for the best possible measuring quality.



Water Resistant

Ability to work in wet canals.

NEO measurement	Approx. 43 x 46 x 16 mm (24g)
File hook length	Approx. 200 mm (w/hook tip)
Lip hook length	Approx. 56 mm
Lip hook diameter	Approx. 2.15 mm
Probe cord length	Approx. 70 cm
Extension cord length	Approx. 60 cm
Stand measurement	Approx. 43 x 67 x 62 mm (52g)

Contact your local **dealer**.

